



AI INTELLIGENT BODY SCAN

A multi-factor setting system and intelligent algorithm are used to provide more comprehensive analysis, more recorded count signals, faster responses and more accurate control when scanning the body and regulating massage techniques.

DOUBLE MECHANISM

Working in combination with the specialised 2D massage robot, the 3D massage robot delivers a synchronised massage that is akin to being massaged with four hands at the same time.

FLEXIBLE TRACK

Because of the long massage path, the flexible track can open into an arch shape to lift the back, which truly achieves Thai stretching massage in a bridge posture.

THE AIR PRESSURE INTENSITY IS ADJUSTABLE

Complete relaxation from the shoulders to the feet with a total of 38 airbags and adjustable air pressure intensity.

DUAL WARMING FUNCTION


Graphene is the newest, strongest, thinnest, most conductive material that is 100% carbon. It surpasses all other known materials as a conductor of heat. The far infrared rays it produces are considered an ideal therapy for the natural prevention of health problems and also promote blood circulation when combined with carbon warmth to the back.

STEREO SPEAKER

Amidst all the attractive features, the armchair also has stereo speakers for deep relaxation of body and soul. The NICAM stereo speakers offer a unique experience of listening pleasure.

LED TOUCHSCREEN


A touch screen featuring the latest UI operating system visualises the massage process. The massage programme can be controlled with just a single button.



BEING IN GOOD CARE AND WELL UNDERSTOOD – YOUR SMART RETREAT



DOUBLE MECHANISM AND A FLEXIBLE TRACK FOR EVEN BETTER COMFORT



AVAILABLE COLOURS



TECHNICAL DATA

Dimensions upright	81,5 x 167 x 124,5 cm
Dimensions reclining	81,5 x 180-195 x 98,5 cm
Weight	159 kg
Voltage	220-240V ~ 50/60Hz
Nominal power	180 W
Rated runtime	30 minutes
braintronics®	21 minutes
Certificates	CE EAC



Width of seating area 46 cm | Width of shoulder area (inside) 55 cm



Casada International GmbH
Obermeiers Feld 3
33104 Paderborn | Germany
www.casada.com

time to relax



TITAN

MASSAGECHAIR



SPECIALS



BRAINTRONICS®

Braintronics® technology synchronises your brainwaves through audio stimulation. Combined with a specially developed massage programme your stress-ridden body will be able to relax. No matter whether you'd like to relax, dream or learn – anything is possible.



SHIATSU
Shiatsu (finger press) is a form of body therapy originating in Japan. The massage technique consists of soft, rhythmic, far-reaching stretches and rotations. The aim of Shiatsu is to stimulate individual parts of the body and to mobilise the muscles.



3D MASSAGE
In this new technique, the massage robot not only moves upwards and downwards, but also forwards and backwards.



2D MASSAGE
The intelligent 2D massage mechanism allows the massage chair to measure your back to determine a more exact position of the sections of the spine.



KNEADING
Skin and muscles are either clasped and kneaded between the thumb and the index finger or with both hands. This massage technique is used in particular to ease tension.



TAPPING
Short, tapping movements are performed with the edge of the hand, the palm or the fist. This promotes blood circulation in the skin and softens tense muscles. If tapping is performed at lung height, this can improve the release of mucus in the lung.



KNEAD & TAP (DUAL ACTION)
The techniques of kneading and tapping are combined in order to relieve tension and stimulate blood circulation at the same time.



ROLLING
Pleasant alternation between relaxation and releasing the muscles – also perfectly suited for a final massage, for relaxation and recovery of the muscles.



CARBON HEATING
Deeply penetrating Carbon infrared heat, blood vessels expand, the blood flow increases and the muscle tone is reduced. In addition, it carries out a highly beneficial effect on nerve cells, which are responsible for reporting pain to the brain.



SOLE ROLLER
The foot massage stimulates the nerve endings and certain organs of the body, which also relaxes the muscles.



S-LINE
With this automatic program the massage heads are working very close to your body, they massage the whole back down to your buttocks. The S-shaped mechanic is constructed after the natural anatomy of the spine.



AIR COMPRESSION
Rising and receding airbags have a pumping motion which produces natural muscle reflexes. This is particularly beneficial in the calf and foot area.



ZERO SPACE
For the lying-position the massage chair moves forward and therefore requires very little distance to the wall behind. Space saving and comfortable!



ZERO GRAVITY
You put yourself in the position of an astronaut – that "weightless feeling".



AUTOMATIC
Specifically developed automatic programmes for the whole body, which offer a big range of different massage techniques. Additionally you can also turn on the heating function, if you'd like.



BACK AREA
An individual applicable to any size massage area.



SHOULDER GRASP
Targeted kneading massage in the shoulder/neck area – ideal for easing tension and pain in this area.



RECLINING POSITION
The adjustability of the backrest and footrest contributes to personal well-being of the body position.



ADJUST SHOULDER
The massage heads can be moved upward or downward to a suitable position in the shoulder area.



AIR INTENSITY
Adjustable airbag massage intensity.



AIR AREA
Whole body airbag massage. Different areas are individually selectable.



SPOT
The massage heads can massage a particular spot.



SPEED
Massage speed is adjustable.



WIDTH
The distance between the both massage heads can be set (narrow, medium and wide).



MUSIC
Further relaxation and fun with the music function.



BLUETOOTH
Wireless data transmission.



HUMAN HANDS FEELING
Specially designed massage heads imitate the hand motion of a professional masseur. The maximum stretch space is 6,5 cm, and the maximum stretch angle is 41 degrees, so the mechanical hands can do effective massage to human cervical vertebra thoracic vertebra and lumbar vertebra.



USB-CHARGER
While you relax, the integrated USB port can supply, for example, your phone with power.



QUICK START
With one touch the quick start buttons in your armrest allow you to control the most important functions.



ARTHROSE TRACTION
Arthrosis Traction consists of gentle stretching pressure. It has a stimulating effect on the cartilages, which enhances the transport of bodily fluids, nutrients absorption and can alleviate symptoms of arthrosis.



WAIST STRETCH
Airbags in the lower area of the back ensure gentle stretching of the pelvis. Effectively counteracts tensions.



L-SHAPE
You can expect a profound full body massage from the neck down to the buttocks.



ADJUSTABLE FOOTREST
The footrest can be extended, this way it is adaptable to different body sizes.



VOLUME
Adjusting the volume.



BACK STRETCH ROLLING
Targeted stretch massages activate the body and have an invigorating effect on the whole body.



ACUPRESSURE POINTS
Key acupressure points are massaged.



REFLEXOLOGY
Reflexology regenerates foot muscles and has a harmonizing effect on almost all organs inside your body.



DUAL SENSOR BODY SCAN
The massage chair performs a body scan to capture the main areas of the neck and back to ensure optimal pressure throughout the massage. The scan technology detects the height of your shoulders, that can be modified after the scan, if required.



SLIM LCD TOUCH SCREEN REMOTE CONTROL
Using the touch screen, you can select your desired massage comfortably and adjust all other settings. You can navigate through the menu by clicking the buttons to select Neck/Shoulder, Lower Back, Auto, Manual, Advance or Air Massage. Optionally, you can set the massage time to a period of 5 minutes, starting from 5 minutes up to 30 minutes.



165° RECLINING POSITION:
Electric, adjustable backrest with an angle of up to 165 degrees.



2 MASSAGEROBOTS
Two massaging mechanisms working separately of each other.



VOICE CONTROL
Use the voice control to activate selected functions.



GRAPHENE WARMTH IN THE LEG AREA
This switches on with the heating function in the lumbar area. The Graphene technology heats the leg area within a few minutes.

braintronics®

brainwave stimulation system

The human brain is a complex organ consisting of millions of little nerve cells, called neurons. Each cell constantly produces tiny electrical impulses. Every thought, every sound and every impression, everything is transmitted and communicated by those small electrical impulses. Electricity is, so to say, the language of the brain. The sum of those impulses can be detected and diagrammed, e.g. by using an EEG. Charted on a graph this electrical activity forms a wave-like pattern, known as brainwaves.

With braintronics®, combining audio stimulation and guided meditations with a soothing massage, these brainwaves can be stimulated in order to reach different states of relaxation in the brain. Whether you want to sustainably enhance your concentration, or just relax and pamper yourself – everything is possible with a massage chair from the Casada family.



CASADA
HEALTHCARE



TITAN